



## Ultra-Pure BlasTaq™ 2X PCR MasterMix

### Cat. No. G885

Store at -20°C.

### Product Description

**Ultra-Pure BlasTaq™ 2X PCR MasterMix** is a ready-to-use MasterMix containing **abm's** Ultra-Pure BlasTaq™ DNA Polymerase in a uniquely-formulated buffer with gel loading dye. This strategically-engineered, next generation Taq Polymerase provides **rapid extension rates and robust performance**. With specialized reaction conditions, this polymerase provides increased processivity, yields, and sensitivity, while shortening reaction times by up to 70%, compared to wild-type Taq DNA polymerase. BlasTaq™ has 5'-3' polymerase and 5'-3' exonuclease activities, lacks 3'-5' exonuclease activity, and produces 3'-dA-tailed amplicons. PCR products made with BlasTaq™ can be used with TA cloning vectors. The enzyme is subject to a rigorous multi-step purification protocol using physical, chemical, and enzymatic methods for maximum removal of contaminating genomic DNA. Quality control involves a PCR DNA Contamination Test with three different primer sets against host DNA.

| Product Component                                 | Quantity          | Part No. |
|---------------------------------------------------|-------------------|----------|
| Ultra-Pure BlasTaq™ 2X PCR MasterMix <sup>1</sup> | 100 rxn (1.25 ml) | G885     |

<sup>1</sup> Buffer contains 1.5 mM Mg<sup>2+</sup>.

### Protocol

1. Mix individual components before use and assemble reaction on ice.

| Product Component                    | Volume                        |
|--------------------------------------|-------------------------------|
| 2X Ultra-Pure BlasTaq™ PCR MasterMix | 12.5 µl                       |
| Forward Primer (10 µM)               | 1 µl                          |
| Reverse Primer (10 µM)               | 1 µl                          |
| Template DNA                         | Variable (100 ng genomic DNA) |
| Nuclease-free H <sub>2</sub> O       | Up to 25 µl                   |

2. Gently mix the reaction components, and briefly centrifuge. Run thermocycling conditions for standard PCR:

| Step                              | Temperature       | Time      |
|-----------------------------------|-------------------|-----------|
| Initial Denaturation <sup>2</sup> | 95°C              | 3 min     |
| 25 – 35 Cycles                    | 95°C              | 15 sec    |
|                                   | 60°C <sup>3</sup> | 15 sec    |
|                                   | 72°C              | 15 sec/kb |
| Final Extension                   | 72°C              | 1 min     |

<sup>2</sup>For most applications, an initial 3 minute denaturation step at 95°C is sufficient. Increase to 5 minutes for high-GC or difficult templates.

<sup>3</sup>BlasTaq™'s PCR buffer allows for primer annealing at 60°C for most primers and adjust only if needed.

3. After PCR, maintain the reaction at 4°C or store at -20°C until use.
4. Analyze the amplification products by agarose gel electrophoresis.
5. Visualize by ethidium bromide or SafeView™ (Cat No. **G108**) staining.

### General Notes

- Specialized buffer for higher yields, sensitivity, and specificity compared to wild-type Taq polymerase.
- Decrease reaction times by 70% using specialized protocol.
- For optimal efficiency, use a 25 µl reaction volume.
- We recommend assembling no template controls in a DNA-free environment to prevent contamination. Nuclease-free pipette tips and PCR tubes will also improve the results.