

contamination risk.

13. The internal threaded screw cap and the silicone O-ring adopt a two-color injection molding process to ensure the sealing of the cap and the vial.

Warning

1. Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquefied gas.
2. Always use appropriate safety equipment when removing vials from cryogenic storage.

Product Range**Cryogenic Vial, Package**

Cat. No.	Product Description	/Pack	/Case
618001	0.5 mL External Thread	50	2,000
606001	1.2 mL External Thread	50	2,000
607001	2.0 mL External Thread	50	2,000
624001	3.0 mL External Thread	50	1,000
608001	4.0 mL External Thread	50	1,000
609001	5.0 mL External Thread	50	1,000
618101	0.5 mL Internal Thread	50	2,000
606101	1.2 mL Internal Thread	50	2,000
607101	2.0 mL Internal Thread	50	2,000
624101	3.0 mL Internal Thread	50	1,000
608101	4.0 mL Internal Thread	50	1,000
609101	5.0 mL Internal Thread	50	1,000

Cryogenic Vial, Racked Package

Cat. No.	Product Description	/Pack	/Case
618102	0.5 mL Internal Thread	100	1600
606102	1.2 mL Internal Thread	100	1400
607102	2.0 mL Internal Thread	100	1200
618003	0.5 mL External Thread	81	1296
606003	1.2 mL External Thread	81	1134
607003	2.0 mL External Thread	81	972
624003	3.0 mL External Thread	81	848

Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi,
Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Overseas


NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab
Emirates)

608003	4.0 mL External Thread	81	648
609003	5.0 mL External Thread	81	486
618103	0.5 mL Internal Thread	81	1296
606103	1.2 mL Internal Thread	81	1134
607103	2.0 mL Internal Thread	81	972
624103	3.0 mL Internal Thread	81	648
608103	4.0 mL Internal Thread	81	648
609103	5.0 mL Internal Thread	81	486

Technical Drawing


External Thread

0.5ml



Internal Thread

0.5ml



Head Office


No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)


External Thread

1.2ml




Internal Thread

1.2ml




External Thread

2ml



Internal Thread

2ml



Head Office


No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi,
 Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab
 Emirates)


External Thread

3ml




Internal Thread

3ml




External Thread

4ml



Internal Thread

4ml



Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)


External Thread

5ml



Internal Thread

5ml



Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)